

# BEFOUR

MEASURING NEW HEIGHTS

1-800-367-7109

www.befour.com

## DISPLAY FEATURES

Measurements are displayed in bright LED digits in Inches, Feet and Inches, or Centimeters.



Inches



Feet and Inches



Centimeters

## Flexible Measurement Range

Measures up to 91 inches (230 cm)



## Model: HTR-101 Wall Mount Digital Stadiometer

The HTR101 Height Rod features fast and accurate digital performance at an attractive price. With its user-selectable 51" measurement range, the versatile HTR101 can measure the height of a small child to that of an NBA player. The smooth glide headpiece features one-button operation, bright LED digits that can be easily seen even at its maximum height, and a unique auto-hold function that retains the measurement display even as it is raised out of the way.

The HTR101 stadiometer has an accuracy of  $\pm 0.1$  in (0.1 cm) and easily converts between three display formats – inches, feet & inches, and centimeters. With no wires to get in the way, the HTR101 is ideal for clinics, hospitals, schools, and sport and fitness centers. Powered by (2) "AAA" alkaline batteries (included), and with its automatic turn off, the HTR101 is good for up to 10,000 measurements.



### FEATURES/SPECIFICATIONS

- Fast and accurate measurement
- View height in inches - feet and inches - centimeters
- Measures up to 91 in (7 ft & 7 in / 230 cm)
- Effortless sliding headpiece
- Flexible measurement range
- Bright LED digits

Befour Inc.

102 Progress Drive  
Saukville, WI 53080

Phone  
800-367-7109

Web  
www.befour.com

Fax  
262-284-5966

E-mail  
sales@befour.com

